HappyCells.Bio

Recombinant Human Interleukin 6 (IL-6) | SKU: RDRCP03

Datasheet:

Product name:	Recombinant Human Interleukin 6 (IL-6) Interleukin-6 (30-212 aa, N-terminal His-T7 tag) Human Protein		
Description:	Recombinant human IL-6 protein fused to His-T7 tag at N-terminus (30-212 aa, His-tag), was expressed in E. coli and purified by using conventional chromatography techniques		
Synonyms:	BSF2, HGF, HSF, IL-6, interferon, beta 2		
UniProt ID:	P05231.1		
Species Expression host Tag:	Human <i>E. coli</i> His-T7 tag		
Expression cDNA Clone or AA Sequence:	MGSSHHHHHHSSGLVPRGSHMASMTGGQQMGRMVPPGEDSKDVAAPHRQPLTSS ERIDKQIRYILDGISALRKETCNKSNMCESSKEALAENNLNLPKMAEKDGCFQSGFNEET CLVKIITGLLEFEVYLEYLQNRFESSEEQARAVMSTKVLIQFLQKKAKNLDAITTPDPTTNA SLLTKLQAQNQWLQDMTTHLILRSFKEFLQSSLRALRQM		
Predicted MW:	~24.3 kDa		
Product Format:	Lyophilized from filtered (0.2 μ m) buffer solution		
Purity:	≥ 95% by SDS–PAGE		
Endotoxin Level:	\leq 1EU/µg as determined by Limulus Amebocyte Lysate (LAL) assay		
Biological Activity:	Measured in a cell proliferation assay using RPMI 8226 cells. The ED50 for this effect is <2 ng/mL corresponding to an activity of >2X10 ⁶ IU/mg		
Storage Temperature:	Store lyophilized at −20°C for up to 1 year		
Stability:	Reconstituted product is stable for 2 weeks at 2–8°C or as aliquots with a carrier protein at –20°C for longer. Avoid multiple freeze-thaw.		
Product details:	R*NR** 100 75 50 37 20 15 10	IL-6	
		Batch no.	DIL6P01/24
		Quantity	10 µg
	*R - Reduced **NR - Non Reduced	Manufacturing	August 2024

RESOLVE BIOTECH PRIVATE LIMITED • 4th Floor, Jet Prime, 39-B, Suren Road, Andheri (East), Mumbai - 400 093 • +91.22.6689.2400 | +91.22.6689.2401 | sales@resolvebiotech.com | +91.22.6689.2400 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.66801 | +91.22.6689.2401 | +91.22.6689.2401 | +91.22.